

# SEQUENTIAL CIRCULATOR 3008

HCPCS CODE E0652

*A Gradient, Sequential, Pneumatic Compression Device for...*

## LYMPHEDEMA VENOUS INSUFFICIENCY WOUND CARE



- *Eight Individual Chambers with a Quick Inflate Cycle!*
- *Chambers 1-2, 3-4, 5-6, 7-8 can be individually adjusted*
- *Four Levels of Calibrated Gradient Pressure*


### INTENDED USE

**Model SC-3008 Sequential Circulator** is a gradient, sequential, pneumatic compression device, intended for the primary or adjunctive treatment of primary or secondary Lymphedema. The device is additionally intended for alternate treatment of chronic venous stasis ulcers and associated venous insufficiency, as well as general treatment of swelling of the extremities. The device is intended for home or hospital use.

*Latest Clinical Studies using BIO COMPRESSION SYSTEMS' SEQUENTIAL CIRCULATORS are available upon request*

### Bio Compression Systems Sequential Circulator Model 3008...

- *Reduces Pain*
- *Mimics the Lymphatic System*
- *Promotes Lymphatic Flow by Moving Fluid in the Proper Physiological Direction (distal to proximal)*
- *Stimulates Circulation to Promote Wound Healing*
- *Delivers Continuous Flow*
- *Bilateral Operation*
- *Provides Effective Prophylaxis for Venous Thrombosis*
- *Lightweight • Portable • Easy to Use • Cost-Effective*



*Arm and Hand Therapy*

### Bio Compression Systems' Biocomfort Garments ...

*Made of high quality 200 Denier Nylon Oxford coated with 3 mils of Polyurethane Available in a variety of sizes and styles including*



*Foot and Leg Therapy*

**BIOCOMFORT GARMENT SIZE CHART ON BACK**

### CONTRAINDICATIONS

Pneumatic compression is contraindicated for the patients with the following: 1) Congestive heart failure; 2) Deep Vein Thrombosis; 3) Inflammatory phlebitis or episodes of pulmonary embolism; 4) Infections in the limb, including cellulitis, without appropriate antibiotic coverage; 5) Presence of lymphangiosarcoma  
WARNING: Federal law restricts this device to sale by or on the order of a physician.



*#1 in Compression Therapy*

120 West Commercial Avenue, Moonachie, NJ 07074

Phone: 800-888-0908 Fax: 201-939-4503

E-mail: [biosystems@biocompression.com](mailto:biosystems@biocompression.com)

Website: [www.biocompression.com](http://www.biocompression.com)

# SPECIFICATIONS

## MODEL 3008

**HCPCS CODE E0652**

Segmental, pneumatic  
compression device  
with calibrated gradient  
pressure and 8  
separate chambers.

Dimensions H x W x D: 5.5 x 12 x 8  
Weight: 7 lbs.  
Pressure Range: 0-125 mm Hg  
Electrical Rating: 120 VAC, 60 Hz, .5A

Cycle Time:  
Inflation: 46 secs.  
Deflation: 6 secs.  
Warranty: 3 years

## BIOCOMFORT GARMENTS SIZE CHART *Now in INCHES/CENTIMETERS*

<u>Upper Extremity</u>	<i>1</i>	<i>2</i>	<i>3</i>	<i>4</i>	<i>5</i>	<i>6</i>	<i>6</i> Axilla to Fingertips Length
	<u>Axilla</u>	<u>Bicep</u>	<u>Elbow</u>	<u>Forearm</u>	<u>Wrist</u>		
<u><i>Small Arm</i></u> G8-3035-S	27"/69.6	26"/66.0	24"/61.0	23"/58.4	19.5"/49.5		26.5"/67.3
<u><i>Medium Arm</i></u> G8-3035-M	27"/69.6	26"/66.0	24"/61.0	23"/58.4	19.5"/49.5		31.0"/78.7
<u><i>Large Arm</i></u> G8-3035-L	27"/69.6	26"/66.0	24"/61.0	23"/58.4	19.5"/49.5		35.5"/91.2
<u><i>Arm w/ Shldr</i></u> G8-3035-SH	27"/69.6	26"/66.0	24"/61.0	23"/58.4	19.5"/49.5		31+8"/78.7+20.3
<u>Lower Extremity</u>	<i>1</i>	<i>2</i>	<i>3</i>	<i>4</i>	<i>5</i>	<i>6</i>	<i>7</i>
	<u>Thigh below Crotch</u>	<u>Mid Thigh</u>	<u>Mid Knee</u>	<u>Mid Calf</u>	<u>Ankle</u>	<u>Arch (Instep)</u>	<u>Crotch to Heel</u>
<u><i>Small Leg</i></u> G8-3045-S	32"/81.3	31"/78.7	25"/63.5	24.5"/62.2	24.5"/62.2	21.5"/54.6	27"/69.6
<u><i>Medium Leg</i></u> G8-3045-M	32"/81.3	31"/78.7	27"/68.6	25"/63.5	24.5"/62.2	21.5"/54.6	32"/81.3
<u><i>Large Leg</i></u> G8-3045-L	32"/81.3	31"/78.7	29"/73.7	25.5/ 64.8	24.5"/62.2	21.5"/ 54.6	36"/91.4
<u><i>Small Leg Adjustable</i></u> <i>(Expanded Position)</i> A8-3045-S	45.5"/115.6	42"/106.7	33.5"/85.1	31"/78.8	24.5"/62.2	21.5"/54.6	27"/69.6
<u><i>Medium Leg Adjustable</i></u> <i>(Expanded Position)</i> A8-3045-M	45.5"/115.6	42"/106.7	35.5"/90.2	31.5"/80	24.5"/62.2	21.5"/54.6	32"/81.3
<u><i>Large Leg Adjustable</i></u> <i>(Expanded Position)</i>	45.5"/115.6	42"/106.7	37.5"/95.3	32"/81.3	24.5"/62.2	21.5"/54.6	36"/91.4

NUMBERS IN ITALICS CORRESPOND TO CUSTOM ORDER FORMS

GARMENTS-SL ISSUE 10-12-05 REV. 10-03-06-F